



## ALS CellCelector™

The ALS CellCelector is a freely configurable tool for the automated transfer of single cells and cell colonies. The patented harvest process (pat. US 9,822,331 B2) allows a gentle cell uptake directly from the culture plate without pre-treatment. Using this technology, highest survival rate and cell integrity are guaranteed.

[Our new developed laminar flow cabinet](#) allows you to work permanently under sterilizable and defined environmental conditions (temperature, CO<sub>2</sub>).

A controlled cell harvest has an increasing importance for the applications of the rapidly expanding life science sector.

